





Linewidth Analyzer

High resolution lineshape spectra and ultra sensitive noise analyzers for narrow and broadband lasers The HighFinesse Linewidth Analyzers of the LWA-100k series are high-end devices for measuring and analyzing highly resolved lineshape and frequency noise spectra of narrow- and broadband lasers. A combination of two interferometric based measurement modes and specifically tailored optics and electronics enable the LWA-100k to quickly acquire various information about the laser source.

## The main features are:

- Large wavelength intervals in the visible and near-infrared regime
- Spectral lineshape and frequency noise analysis
- Intrinsic linewidth measurement range from 1 kHz up to 1 MHz
- Effective linewidth measurement range from 100 kHz up to 300 MHz
- Frequency noise sensitivity below 100 Hz/√Hz for frequencies above 100 kHz
- High sample rate up to 30 MSa/s and fast acquisition rates up to real-time
- Linewidth accuracy down to 20 kHz
- No reference source required

The LWA-100k is perfectly suited for laser development and adjustment. In combination with a HighFinesse wavemeter high resolution classification of lasers and laser systems can be performed.



The HighFinesse instruments of the LWA-1k series are the ultimate high-end devices for measuring, analyzing and controlling frequency and intensity noise of lasers. The superb sensitivity of the LWA-1k series is achieved by combining an interferometric working principle with high-end optical and electronic components..

## The main features are:


- Frequency noise, spectral lineshape and relative intensity noise (RIN) analysis with evaluation of intrinsic and combined intrinsic and effective linewidth
- Intrinsic linewidth measurement range down to 350 Hz
- Effective linewidth measurement range down to 2 kHz
- Frequency noise spectrum sensitivity down to 10 Hz/√ Hz with a dynamic range of 50 dB between 10 Hz and 10 MHz
- RIN down to -150 dB/Hz
- Extremely robust against acoustic noise
- Error signal generator for further linewidth, frequency noise or RIN reduction
- Powerful tool for a detailed analysis of noise sources like servo bumps, frequency drifts, power supply noise and acoustics.





LWA-100k Series LWA-1k-1550

The LWA-100k analyzes both, very narrow laser lines down to 100 kHz as well as broader bands up to 300 MHz. This is achieved by two distinct measurement modes. In the broadband-mode, a frequency-domain intensity spectrum is evaluated by determining the full width at half maximum (FWHM) of a curve fitted to the experimental data. In the narrowband-mode, time-domain recorded laser frequency deviations are analyzed yielding both a single-sided noise spectrum and a two-sided lineshape spectrum. Once connected to the PC via ethernet, the LWA-100k is controlled

by an intuitively usable software interface that automatically evaluates and presents the data to the user. By monitoring the changes in the lineshape parameters over time, the stability of the light source can be determined. The LWA is ideal for optimizing the stability of laser setups by revealing mechanical, acoustic and electronic noise. Also, a high spectral resolution of 100 kHz in the broadband-mode allows supervising the single-mode operation of the laser and a free spectral range of 1 GHz enables the identification of sidebands.








The LWA-1k featuring an improved resolution enables the analyzation of laser with very narrow linewidths down to 2 kHz. The analyzer unit automatically generates an output signal that is evaluated by the digitizer unit, connected to a PC via USB. An intuitively usable software interface performing the analysis provides the frequency noise related plots and values, such as the frequency noise spectrum and the lineshape spectrum together with the results of intrinsic and combined linewidth.

The frequency noise spectrum is a powerful tool for a detailed analysis of noise sources, servo bumps or frequency drifts enabling e.g. further laser locking improvements.

Additionally, the output signal provided by the analyzer unit of the LWA-1k series can be used as an error signal for further linewidth or RIN reduction in combination with an appropriate hardware such as a PID-controller.

As an example, the graphs on the lefts show frequency noise reduction of an ECDL with 1550 nm in the lab.



Spectrogram of beat between an ECDL and a fiber laser with (black) and without (red) external PID-lock of the ECDL to the LWA-1k confirming the lineshape spectrum in the graph to left.

#### Technical Data

LWA-100k series	Unit	400	500	750	980	1550			
Wavelength range		380 – 430	430 – 660	615 - 885	825 – 1200	1200 - 1700	1) Best performance with		
Required input power <sup>1)</sup>	mW			2 - 5			maximum input power. Noise sensitivity scales inversly with input power.		
Analyzer Unit							Effective linewidth:     Combination of intrinsic     linewidth and additional		
Laser type		CW and single-mode				broadening mechanisms (thermal, electronical and acoustic noise). Determina-			
Input type		FC/APC fiber							
Broadband mode specification							tion by β-separation method (noise density spectrum)		
Optical frequency resolution	Hz	2 M				or curvefitting procedure (lineshape spectrum). 3) Intrinsic linewidth: Limited by fundamental			
Free spectral range	Hz	1 G							
Effective linewidth <sup>2)</sup> range (FWHM)	Hz	2 M – 300 M							
Spectral and frequency noise specification	ns						quantum processes and laser design. Determined by the		
Optical frequency resolution	Hz					noise floor (white noise) of the frequency noise spectrum and calculated by: noise density (in Hz²/Hz) times π (rule of thumb). This value is most commonly denoted as "laser linewidth" by laser			
Frequency noise bandwidth	Hz	100 – 1 M							
Frequency noise sensitivity	Hz/√Hz	<100 (@ > 100 kHz)							
Intrinsic linewidth range <sup>3)</sup>	Hz	1 k - 1 M							
Effective linewidth range (β-separation) 2)	Hz			100 k – 1 M			manufacturer.		
Lineshape spectrum									
Effective linewidth 2) range (FWHM)	Hz			100 k – 1 M			-		
Optical frequency resolution	Hz	100 k							
Miscellaneous									
Dimensions	mm			150 × 280 × 7	9		-		
Weight	kg	5							
Digitizer module									
Sample rate	Sa/s			31 M (max.)			-		
Resolution	bits			16			-		
Acquisition time	S			0.1 (typ.)			-		
Evaluation time	S			<1 (typ.)			-		
Miscellaneous									
Communication				Ethernet			-		
Dimensions	mm		;	357 × 112 × 14	15		-		
Weight	kg			8			-		
Software									
Operating system			Microsoft®	Windows® (7	– 10), 64 Bit		-		
CPU (minimum)			Intel® C	ore™ i5 or eq	uivalent		_		
Memory (minimum)				8 GB			_		
Ports				Ethernet			-		
Graphical Evaluation options					noise density distribution (				

#### Technical Data

**Graphical Evaluation options** 

LWA-1k series	Unit	LWA-1k 780	LWA-1k 1550	
Wavelength range	nm	780 – 1020	1530 - 1565	Best performance with     maximum input power.
Required input power <sup>1)</sup>	mW	1- 10	0.5 - 5	Noise sensitivity scales
Input power stability	<u></u> %	±5		inversly with input power.
Analyzer Unit				Effective linewidth:     Combination of intrinsic     linewidth and additional     broadening mechanisms     (thermal, electronical and
Laser type		CW and single-mode		acoustic noise). Determina-
Input type		PM-FC/APC fiber		tion by β-separation method (noise density spectrum)
Spectral and frequency noise specifications				or curvefitting procedure (lineshape spectrum).
Optical frequency resolution	Hz	8 k	1 k	_
Frequency noise bandwidth	Hz	10 – 10 M		3) Intrinsic linewidth: Limited
Frequency noise sensitivity	Hz/√Hz	< 50 – 10 M	< 10 - 10 M	by fundamental quantum processes and laser design.
Intrinsic linewidth range 3)	Hz	< 8 k	<350	Determined by the noise floor (white noise) of the
Effective linewidth range (β-separation) 2)	Hz	< 20 k - 10 M	< 2 k – 10 M	_ frequency noise spectrum
Relative intensity noise limit	dB/Hz		150	and calculated by: noise  density (in Hz²/Hz) times π
Lineshape specifications				(rule of thumb). This value is most commonly denoted
Effective linewidth 2) range (FWHM)	Hz	< 20 k - 4 M	< 2 k - 4 M	as "laser linewidth" by laser
Optical frequency resolution	Hz	20 k	2 k	manufacturer.
Dynamic range	dB	!	50	<ul> <li>4) Linewidth reduction/control:</li> <li>Analog output as error signal</li> </ul>
Miscellaneous				for use in combination with  PID controller (not included)
Interface		USB 2.0 Type B		for frequency noise or RIN
Analog Output/error signal <sup>4)</sup>		BNC $\pm$ 7.5 (50 $\Omega$ ) $\pm$ 15 (high impedance) V, single ended		reduction.
Cutoff (highpass filter)	Hz	10, 1 k, 10 k, 100 k		_
Dimensions	mm	220 × 334 × 96		_
Weight	kg	8		_
Digitizer module				
Sample rate	Sa/s	62.5 (max.)		_
Resolution	bits	16		_
Acquisition time	S	0.1 (typ.)		_
Evaluation time	S	<1 (typ.)		_
Miscellaneous				
Communication		USB 3.0 Type B		_
Dimensions	mm	210 × 200 × 74		_
Weight	kg	2		-
Software				
Operating system		Microsoft® Windo	_	
CPU (minimum)		Intel® Core™ i	Linewidth Analyzer · 07-2019  This document provides general	
Memory (minimum)		8	information only and may be	

Frequency noise density graph, lineshape graph,

frequency deviation distribution (histogram)

subject to change at any time

without prior notice.





# Wavelength Meter

HighFinesse/Ångstrom offers sensitive and compact wavelength meters with a large spectral range for high speed measurement of lasers. The optical unit consists of temperature-controlled Fizeaubased interferometers that are read out by photodiode arrays. The high absolute accuracy is achieved by use of solid state, non-moving optics. The optical unit and associated electronics are housed in a compact, thermal casing. The connection to a computer or notebook is realized via a highspeed USB 2.0 port, which allows a high data read-out rate. The analyzing software displays all the interferometer information.



# Spectrometer OSA

HighFinesse/Ångstrom optical spectrometers LSA and HDSA are designed to analyze the multiline or broadband spectrum of light sources like cw and pulsed lasers, gas discharge lamps, super luminescence diodes, semiconductor laser diodes and LEDs. They are suitable to analyze the spectrum of telecom signals, resolve Fabry-Perot modes of a gain chip, and produce a spectral measurement of gas absorption.



# Precision Current Sources

HighFinesse Precision Current Sources have been developed for experiments and quantum technologies in the areas of Cold atom physics and solidstate-physics. The linearly regulated BCS (Bipolar Current Source) and UCS (Unipolar Current Source) series deliver highly stable, low noise source currents for high precision magnetic field control. The current output is floating or is on a user defined potential. Ultrafast response to control signals and trigger functions, clear grounding, connection and signal isolation schemes make the integration of the current sources into complex experimental systems easy.



HighFinesse GmbH Wöhrdstraße 4 72072 Tübingen, Germany







T + 49 (0) 7071 - 53 918 0 F + 49 (0) 7071 - 53 918 99 M info@highfinesse.com





Additional information and distributors: www.highfinesse.com