



Multimode switches

Quick Coupling Specs



HighFinesse
Laser and Electronic Systems

Quick coupling accuracy (with multimode fiber)

WS5	3000 MHz
WS6-600	600 MHz
WS6-200	600 MHz (200 MHz for WS6-200 IR-III)
WS7-60	150 MHz
WS7-30	100 MHz
WS8-30	100 MHz
WS8-10	100 MHz
WS8-2	100 MHz

Further Information

For further technical information please contact:

Pascal Notz

service@highfinesse.de



HighFinesse GmbH
Wöhrdstraße 4
72072 Tübingen, Germany



T + 49 (0) 7071 - 53 918 0
F + 49 (0) 7071 - 53 918 99
M info@highfinesse.com



Additional information
and distributors:
www.highfinesse.com